

Eukaryotic RNA Sequencing

1. Sample Requirements

Library Type	Sample Type	Amount	Volume	Concentration	RNA Integrity Number (Agilent 2100™)	Purity (NanoDrop™)
Eukaryotic RNA-Seq (cDNA library)	Total RNA	≥ 200 ng	≥ 10 μL	≥ 20 ng/μL	≥ 4.0, with smooth base line	OD260/280 ≥ 2.0; OD260/230 ≥ 2.0; No degradation, No contamination
	Total RNA (Blood)	≥ 400 ng	≥ 20 μL	≥ 20 μL	≥ 5.8, with smooth base line	
	Total RNA (Single Cell)	≥ 100 ng	≥ 20 μL	≥ 10 ng/μL	≥ 5.8, with smooth base line	OD260/280 ≥ 2.0 OD260/230 ≥ 2.0; No degradation, no contamination
	Amplified cDNA (double-stranded)	≥ 100 ng	≥ 10 μL	≥ 10 ng/μL	Fragments should be distributed between 400bp-5000bp with main peak at ~2000 bp	
Eukaryotic RNA-Seq (strand specific library)	Total RNA	≥ 400 ng	≥ 20 μL	≥ 20 ng/μL	≥ 5.8, with smooth base line, with smooth base line	OD260/280 ≥ 2.0; OD260/230 ≥ 2.0; No degradation, no contamination

For total RNA less than 100 ng, please contact us for ultra-low input solutions.

2. Sequencing Parameters

Platform	Illumina NovaSeq 6000
Read length	Paired-end 150
Recommended sequencing depth	≥ 20 million read pair per sample for species with reference genome; ≥ 50 million read pairs per sample for species without reference genome (de novo transcriptome assembly projects)
Data quality	Guaranteed ≥ 85% bases with Q30 or higher
Turnaround time	Within 2~3 working weeks from library construction verification to data releasing without bioinformatic analysis. (depending on the sample size);

3. Data Analysis Contents

3.Data Analysis Contents

Standard analysis
Data filtering
Transcriptome assembly & Gene functional annotation (only for species without reference genome)
Mapping to reference genome/assembled genome
Gene expression quantification & Differential expressed genes profiling & Enrichment analysis
Protein-Protein Interaction (PPI) analysis
Transcription factors functional annotation analysis
Oncogene functional annotation analysis
SNP & InDel analysis
Alternative splicing analysis
Fusion gene prediction (Only for tumor sample and cancer cell line)